

Flexible Video CHOLEDOCHOSCOPE CH-1349



CH-1349

CMOS + LED Video Choledochoscope



Chip-on-Tip
Video Technology




Bending Up/Down
160°/130°

Photograph and
Video Recording



Compatible with
displays of
mainstream brands



Parameter

Model: CH-1349

Insertion Tube(mm): 4.8

Channel(mm): 2.0

Bending(Up/Down): 160°/130°

Working Length(mm): 380

Field of View: 90°

Image Processor: HV-3101/CV-1080

