

# Flexible Video URETERORENOSCOPE

## UR-1328




# UR-1328

CMOS + LED Video Ureterorenoscope



Chip-on-Tip  
Video Technology



275° bending  
to lower calyx

Left & right 120°  
insertion rotation  
facilitates therapy



Photograph and  
Video Recording

## Parameter

Model: UR-1328

Insertion Tube: 8.5Fr

Bending(Up/Down): 275°/275°

Working Length(mm): 670

Rotation(Left/Right): 120°/120°

Field of View: 120°

