




- Cannula and stylet designed to allow an easier penetration of the bone cortex
- Adjustable depth stop for sternal or paediatric application, which can be removed for iliac crest use
- Luer slip and luer-lock syringes fit securely



Side distal holes VVIE1xxyy/xx(x)

MYELO SLIM Qty/Box 10pcs

VVIE1xxyy/xx(x) for bone marrow harvesting needle

CODE	SIZE
VVI11410/50 VVIE11410/50	14Gx30mm
VVI11430/70 VVIE11430/70	14Gx50mm
VVI11510/50 VVIE11510/50	15Gx30mm
VVI11530/70 VVIE11530/70	15Gx50mm
VVI11610/50 VVIE11610/50	16Gx30mm
VVI11630/70 VVIE11630/70	16Gx50mm
VVI11810/50 VVIE11810/50	18Gx30mm

On request Vigeo is available to provide needles with different diameter and length