

AsPEKT A703/712 digital recorders with LCD screen

District Distring District District <td

AsPEKT A703/712 recorders:

- presentation of the ECG signal at the LCD screen,
- easily detachable patient cable,
- ▶ SD card support,
- Bluetooth connection with PC.

AsPEKT A703/712 digital recorders



AsPEKT A703/712 recorders:

- little portable dimensions,
- recording of 3 channels (AsPEKT 703)/ 12 channels (AsPEKT 712),
- SD card support,
- Bluetooth connection with PC.

Disclaimer: all materials, descriptions and images included in this brochure are for informative purposes only. Product images may differ from actual product appearance.

Aspel SA

os. Henryka Sienkiewicza 33, 32-080 Zabierzów, Poland phone: +48 12 285 22 22; fax: +48 12 285 30 30 mail to: export@aspel.com.pl www.aspel.eu



IN ISO 9001 IN ISO 13485 MDD 93/42/EEC Free Sale Certificates



10.2017

ASPEL INNOVATION EXPERIENCE QUALITY

HolCARD 24W

HolCARD 24W is an easy to operate software enabling efficient moving between ECG registration and analysis. In cooperation with AsPEKT recorders it creates a system to be installed on an user's PC and is dedicated to carry out both quick, basis and detailed, advanced analysis of ECG Holter examination.

🖶 HolCARD 24W (c) Aspel - Patient: Adam Nowak Examination: 13.06.10 09:00	
Examination Analysis Review Report Settings Window Help Set lock Tem 12 TE TE TE C Set UP Lock Review Tem 20 C Set UP Lock Review Tem 20 C Set UP Lock Review Report Settings Window Help	Artef: 0.00%
	11111111111111
	p-p-d-d-d-d-d-d-d-d-d-
Martin and the second of the s	~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~
	-barded and have been a start of the second st
06:57 - A.A.A.A.A.A.A.A.A.A.A.A.A.A.A.A.A.A.A	+++++++++++++
	-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h
	-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h-h
	hhhhhhhhhh
	-hhhhhhhhhhhh
~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~	had
How hat	- Add to the hard
- A-A-A-A-A-A-A-A-A-A-A-A-A-A-A-A-A-A-A	n min
855 K ← → ≫ 8 12 20 ± 2501 30: 452 60: 🔛 → 🖨 Chan 3 •	
06:54:00	
man and a start and a start	
	~-
3 The hard and a state of the s	
RR HR K ← → >> 5. 12. 22	
ASPEL - System Hother EKG	NUM
	and the second

- estimation of arrhythmias and of ST segment
- analyses of ECG registration in prospective and retrospective modes
- detection, reviewing and compiling quantitative reports within the scope of: VE, single·SVE, PSVT and many others
- results' presentation in a trend format, likewise ST trends in 3-D option
- > analyses of arrhythmias of ventricular and supraventricular origins in registrations till 48 hours
- easy edition of all ventricular and supraventricular impulses owing to the used class editor function
- independent analyses of ST segment of each ECG channel for chosen by a user ST points of reference; J point as well as R-R base line

Disclaimer: all materials, descriptions and images included in this brochure are for informative purposes only. Product images may differ from actual product appearance.

Aspel SA



www.aspel.eu







10.2017