



ASCO
SINCE : 1948

Apothecaries' Sundries Manufacturing Co.

CATALOGUE # LED15


www.ascomedical.com



MEMORY FUNCTION



ENDO MODE



LIGHT INTENSITY

80,000 + 1,40,000 Lux $\pm 10\%$

COLOUR TEMPERATURE



3500-5000 K (ADJUSTABLE)

NUMBER OF LED



57 + 95

TECHNICAL FEATURES

Size of Light Field	12-30 cm
Colour Reduction Index	93 RA
LED Life	>50000 Hrs.
Diameter of Light	500mm x 2
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Focusing	Adjustable
Battery Backup	Optional

SPECIAL FEATURES

Memory Function, High quality sterilizable handle, With protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

M MEMORY FUNCTION

E ENDO MODE



LIGHT INTENSITY



80,000 + 80,000 Lux \pm 10%

COLOUR TEMPERATURE



3500-5000 K (ADJUSTABLE)

NUMBER OF LED



57 + 57

TECHNICAL FEATURES

Size of Light Field	12-30 cm
Colour Reduction Index	93 RA
LED Life	>50000 Hrs.
Diameter of Light	500mm x 2
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Focusing	Adjustable
Battery Backup	Optional

SPECIAL FEATURES


Memory Function, High quality sterilizable handle, With protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

CSL5.1001

SERIES 5

M MEMORY FUNCTION

E ENDO MODE



LIGHT INTENSITY



1,40,000 Lux $\pm 10\%$

COLOUR TEMPERATURE



3500-5000 K (ADJUSTABLE)

NUMBER OF LED



95

TECHNICAL FEATURES

Size of Light Field	12-30 cm
Colour Reduction Index	93 RA
LED Life	>50000 Hrs.
Diameter of Light	500mm
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Focusing	Adjustable
Battery Backup	Optional

SPECIAL FEATURES


Memory Function, High quality sterilizable handle, With protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

CSL5.1003

SERIES 5

M MEMORY FUNCTION

E ENDO MODE



LIGHT INTENSITY



80,000 Lux $\pm 10\%$

COLOUR TEMPERATURE



3500-5000 K (ADJUSTABLE)

NUMBER OF LED



57

TECHNICAL FEATURES

Size of Light Field	12-30 cm
Colour Reduction Index	93 RA
LED Life	>50000 Hrs.
Diameter of Light	500mm
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Focusing	Adjustable
Battery Backup	Optional

SPECIAL FEATURES

Memory Function, High quality sterilizable handle, With protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.