



**LED**  
OPERATING  
LIGHTS

**ASCO**  
SINCE : 1948

*Apothecaries' Sundries Manufacturing Co.*

CATALOGUE # LED15

[www.ascomedical.com](http://www.ascomedical.com)



MEMORY FUNCTION



ENDO MODE

**LIGHT INTENSITY**80,000 Lux  $\pm 10\%$ **COLOUR TEMPERATURE**

3500-5000 K (ADJUSTABLE)

**NUMBER OF LED**

57

**TECHNICAL FEATURES**

Size of Light Field	12-30 cm
Colour Reduction Index	93 RA
LED Life	>50000 Hrs.
Diameter of Light	500mm
Brightness Control	Capacitive Touch Panel
Power Supply	220V/50Hz AC
Focusing	Adjustable
Battery Backup	Optional

**➤ SPECIAL FEATURES**

Memory Function, High quality sterilizable handle, With protective fire resistant ESG safety glass for scratch proof, fast disinfection process, protects UV and IR rays and improve light penetration. With special European design for ventilation of heat. Body is made of high quality fire resistant polymer.

**ML5.1003**

Mobile Surgical Light

**ML5.1003B**Mobile Surgical Light  
With Battery Backup