



Description

Hospital Electric Surgical Operation Bed for BT667E

Size: 2100*500*750*1000

- Material: 304stainless steel
- Reverse Trendelenburg:25°
- Lateral Adjustment:20°
- Longitudinal displacement:300mm
- Headrest adjust:- 90°~+45°
- Backrest adjust:- 10°~+60°
- Waist bridge: 100mm
- Leg plate(s): 90°
- Swiveling of leg plate:90°



This operation table is used for the surgeries of chest ,abdomen

,ophthalmology,obstetrics,and urology,etc.It has the following advantages:

- The board can shift longitudinally. It is made of Transparent material material.it can be X-rayed or scan when using along with C-ARM .
- The board is controlled by electric handspike to move up and down, lean from left to right and move back and forth.