



## Description

### 5-Function ICU Electric Hospital Bed for BT605EPZ

Size: 2140\*1000\*440~770mm

Packing size: 223\*107\*50cm(honeycomb box)/1pc

Max Loading: 250KG(550lb) in movement

Function: Backrest tilting From 0° to 75° ±5°

Footrest tilting From 0° to 40° ±5°

Reverse: 14° ±1°

Trendelenburg: 14° ±1°

Sit position function

Hi-lo function

Abdomen auto-extension decompression

Backrest angle indicator

Whole bed angle indicator

Bumper wheels

Handset

#### Standard Accessories:

4pcs 5" Taiwan castors with central brake system,

1 pair brake pedal,

4pcs PP side rails,

2pcs ABS head and foot board,

4pcs urine hooks,

1 set MOTECK motor system,

4pcs bumper wheels,

4pcs angle indicators,

1pc electric cable holder,

4pcs Dust cover

Option: I.V pole, mattress, over bed table, bedside locker