



Description

Hospital Hydraulic Operation Bed for BT668

Size : 2100*500*800~1000

Material : stainless steel

Head function:-90°~ +45°

Back function:-10°~+60 °

Waist function: 100°

Leg function : 90°

Leg abduction:90°

Features:

This table is designed for surgical procedures of the head,neck ,thorax and abdomen ,and extremities as well as for gynaecological ,oto-rhino-laryngological and orthopaedics operations ,

it can be raised or lowered hydraulically by means of an oil pump

various positions can be controlled on the head

it is made of double decker plastic ,X-ray examination can be available during the operation

