



Medical Head Lamp

BT - 410F

BT - 410A

 **bistos**

Medical Head Lamp



BT-410 F

Illumination Uniformity

Multi-lens structure ensure high level of uniform luminance in its light spot



Coaxial Design

Light source is located between two eyes is able to cast least shadow on the surface. Also, the luminance angle moves freely with pivot joint structure.



Interchangeable Battery

Rechargeable battery is interchangeable. With additional battery, no need to worry anymore about remaining power.



Medical Head Lamp



BT-410 A

Adjustable Light Spot

Light spot is adjustable from 30mm to 70mm at distance 300mm.



Astral Lamp

Additional astral Lamp increase more illumination and decrease the shadow on the surface



Loupe

Single lens loupe is light weight and simple but provides 2x magnification.



BISTOS Technical Specification

* Option

Model		BT - 410	
Order code		410F	410A
Illumination	Standard	15,000 - 30,000 lx	30,000 lx
	Astral lamp		50,000 lx
Adjustable Illumination		O	X
Illumination Uniformity		High	Low
Adjustable Light Spot		X	O
Light Spot Size (at 300mm)		70mm	30 - 70mm
Astral Lamp *		X	O
Color Temperature		6,000 K	
LED Life Time		Appox. 50,000 hour	
Loupe *	Type	Single Lens	
	Magnification	2 x	
Battery	Type	Lithium-ion	
	Capacity	3.7V / 2,200 mAh	
	Run time	4 hour	
	Charging time	4 hour	
Weight		190 g	
Band Size		540 - 640mm	
Warranty		1 year	

BIO SIGNAL TOTAL SOLUTION

Bistos Co., Ltd. (Headquarter)

7th Fl., A Bldg., Woolim Lions Valley 5-cha, 302,
Galmachi-ro, Jungwon-gu, Seongnam-si,
Gyeonggi-do, Korea
Tel. +82 (31) 750 0340
Fax. +82 (31) 750 0344

bistos@bistos.co.kr
www.bistos.co.kr

Bistos America Inc. (US Branch)

22941, Triton Way, Suite #242, Laguna Hills, CA., USA 92653
Tel. +1-949 614 8745
Fax. +1-949 614 8745
Site. www.bistosamerica.com
E-mail. info@bistosamericainc.com