



ESSENTIAL LARYNGOSCOPE SETS

Timesco[™]

LARYNGOSCOPE SETS



As one of the largest manufacturers of quality laryngoscope systems in the world we offer a wide selection of pre-configured sets. Sets are available in paediatric, adult and combination versions, in both reusable and single use variants.

Fibre Optic sets are marked in green ● Standard light sets are marked in black ●

All reusable laryngoscope systems are covered by our unrivalled Lifetime Guarantee against faulty materials or manufacture. Consumable parts, such as bulbs and fibre optic cores, are not covered by the lifetime guarantee as they will deteriorate over a period of time.

If you need any further assistance please email export@timesco.com, call +44 (0)1268 297 730 or fax +44 (0)1268 297 803.

OPTIMA CLX



The new fibre optic core in Optima CLX 'clicks' in and out of the block for easier disassembly and assembly. No tools are required to remove it making it a quick and easy process.

The fibre optic core can be replaced at the end of its working life and the blade can be used continuously. The lifetime guarantee on the blade provides a high quality blade that can be used for years.

Supplied with our powerful Optima XL Xenon light handle.



Optima CLX block with removable core

Optima CLX	Handle		Blade No.0		Blade No.1		Blade No.2		Blade No.3		Blade No.4		Bulb	Case
	Small	Medium	Mac	Miller	Mac	Miller	Mac	Miller	Mac	Miller	Mac	Miller		
Standard Set 4092.975.80		•					•		•		•		•	•
Paediatric Set 4092.985.70	•			•		•							•	•
Combination Set 4092.989.95	•	•		•		•	•		•		•		•	•

OPTIMA



Optima fibre optic blades have been our most popular fibre system for almost 20 years, the CLX system above will replace Optima in 2009 and we will slowly be phasing all customers over to CLX. We will continue to supply Optima while stocks are available. Optima has a replaceable fibre optic core which is secured by a grub screw.

All Optima sets are supplied with our powerful Optima Xenon Handles.



Optima Macintosh blade rear profile

Optima	Handle		Blade No.0		Blade No.1		Blade No.2		Blade No.3		Blade No.4		Bulb	Case
	Small	Medium	Mac	Miller	Mac	Miller	Mac	Miller	Mac	Miller	Mac	Miller		
Standard Set 4090.975.80		•					•		•		•		•	•
Paediatric Set 4090.985.70	•			•		•							•	•
Combination Set 4090.989.95	•	•		•		•	•		•		•		•	•

SIRIUS



The Sirius laryngoscope system has been designed with an enclosed fibre optic core. As the fibre bundle is sealed within the blade there are no crevices to trap dirt or body fluids, allowing the blade to be easily and simply decontaminated and cleaned without the need for any disassembly.

Supplied with our powerful Sirius Xenon light handle to provide bright, white illumination. This is then directed upwards into the fibre optic bundle in the blade.



Sirius handle and Macintosh blade.

Sirius	Handle		Blade No.0		Blade No.1		Blade No.2		Blade No.3		Blade No.4		Bulb	Case
	Small	Medium	Mac	Miller	Mac	Miller	Mac	Miller	Mac	Miller	Mac	Miller		
Standard Set 2955.975.80		•					•		•		•		•	•
Paediatric Set 2955.985.70	•			•		•							•	•
Combination Set 2955.989.95	•	•		•		•	•		•		•		•	•

To purchase additional spare bulbs please contact your sales representative.

ORION



Orion is a high quality standard light or conventional Light system. The bulb is housed in the blade and provides bright, focused illumination to the intubation site. Orion blades and handles are fully autoclaveable.

The Orion range is colour coded with black dots on the blades and black bands on the handles. This is for ease of identification that it is a standard light system.



Orion Macintosh blade rear profile

Orion	Handle		Blade No.0		Blade No.1		Blade No.2		Blade No.3		Blade No.4		Bulb	Case
	Small	Medium	Mac	Miller	Mac	Miller	Mac	Miller	Mac	Miller	Mac	Miller		
Standard Set 2905.975.80		•					•		•		•		•	•
Paediatric Set 2905.985.70	•			•		•							•	•
Combination Set 2905.989.95	•	•		•		•	•		•		•		•	•

CALLISTO

2



'No touch' system
with an extended
block lip



Callisto blades are fibre optic single use metal/plastic hybrid blades which are almost identical to reusable blades in feel and performance but at a low cost, disposable price. All blades are compatible with ISO fitting fibre optic handles and are supplied in peel pouches.

The Callisto range of blades has an extended block lip that prevents the tip of the blade coming into contact with the barrel of the handle, preventing cross contamination.

Our Callisto sets are available with either a reusable Optima XL handle or a single use Callisto S plastic handle. Please select the correct handle set when ordering.



Callisto S handle

Callisto	Handle		Blade No.0		Blade No.1		Blade No.2		Blade No.3		Blade No.4		Pouch
	Optima XL	Callisto S	Mac	Miller	Mac	Miller	Mac	Miller	Mac	Miller	Mac	Miller	
Standard Set DS.3940.920.10	•						•		•		•		•
Major Set DS.3940.920.20	•			•		•	•		•		•		•
Standard Set DS.3940.930.10		•					•		•		•		•
Major Set DS.3940.930.20		•		•		•	•		•		•		•

EUROPA

2



Europa Macintosh
blade rear profile



Europa blades are conventional light single use metal/plastic hybrid blades which are almost identical to reusable blades in feel and performance but at a low cost, disposable price. All blades are compatible with ISO fitting conventional light handles and are supplied in peel pouches.

Our Europa sets are available with either a reusable Lunar handle (which is our low cost reusable conventional light handle) or a single use Europa S plastic handle. Please select the correct handle set when ordering.



Europa S handle

Europa	Handle		Blade No.0		Blade No.1		Blade No.2		Blade No.3		Blade No.4		Pouch
	Lunar	Europa S	Mac	Miller	Mac	Miller	Mac	Miller	Mac	Miller	Mac	Miller	
Standard Set DS.2940.920.10	•						•		•		•		•
Major Set DS.2940.920.20	•			•		•	•		•		•		•
Standard Set DS.2940.930.10		•					•		•		•		•
Major Set DS.2940.930.20		•		•		•	•		•		•		•

Timesco

Timesco of London Ltd. Timesco House, 3 Carnival Park, Carnival Close, Basildon, Essex SS14 3WN

T: +44 (0) 1268 297 730 F: +44 (0) 1268 297 803 E: export@timesco.com

www.timesco.com