## autopress / autopress.e

# AUTOMATIC PRESSURE INFUSOR PRECISE PRESSURE SETTING. PRESET PRESSURE.

BIEGLER AUTOPRESS is another component in our range of warming equipment. It can be used either directly with our blood and infusion warmers or as a stand-alone solution where fluids are to be dispensed while maintaining constant pressure in the cuffs.





# autopress / autopress>e

### **AUTOMATIC PRESSURE INFUSOR**



- Preset pressure is maintained in the cuffs .
- Precise pressure setting due to electronic control .
- For use with pressure cuffs 2 x 500, 2 x 1000 or 2 x 3000 ml .
- Mounting on normed railing and infusion stand .
- Compact design .
- Pressure limited to 300 mmHg •
- autopress: Mains operated . autopress • e: Mains operated with up to 3h battery backup
- Pressure cuffs can be emptied quickly and easily •
- Significant reduction of set up time •

### **TECHNICAL DATA**

Device:	a
Type designation:	a
Operating voltage:	
Power consumption:	
Supply type:	
Degree of protection against electric shock:	
Classification (93/42/EEC):	
Operation mode:	
Dimension:	
Weight:	

Accessories:

utomatic pressure infusor

utopress

autopress ► e

100 - 240 V AC, 50-60 Hz

max. 36 VA

mains operated

mains and battery operated

Type CF defibrillationproof

IIa according to Rule 9

continuous operation

100 x 230 x 180 mm



Pressure cuff 500 ml

Pressure cuff 1000 ml Pressure cuff 3000 ml





### **BIEGLER GmbH**

Allhangstrasse 18a, 3001 Mauerbach, Austria Tel: +43 1 979 21 05, Fax: +43 1 979 21 05 16 office@biegler.com, www.biegler.com

