





	Gastrointestinal videoscope		Colonovideoscope	
	OEV-G166	OEV-GP166	OEV-C166	OEV-CT166
Field of View		120° Forward viewing		120° Forward viewing
Depth of View Field		5mm-100mm Fixed Focus		5mm-100mm Fixed Focus
Range of Tip Bending		Up 210°, Down 100° Left 100°, Right 100°		Up 180°, Down 150° Left 150°, Right 150°
Outer Diameter of Distal End	ø 9.8mm	ø9mm	ø 13.2mm	ø 14.3mm
Outer Diameter of Insertion Tube	ø 9.8mm	ø9mm	ø 12.8mm	ø 13.7mm
Inner Diameter of Instrument Channel	ø 2.8mm	ø 2.2mm	ø 3.8mm	ø 4.2mm
Working Length	1025mm	1025mm	1300mm/1500mm/1680mm	
Waterproof	•	•	•	•