

IDS-1000



Standard Composition

IDS-1000

- IDS-1000 Dermatoscope
- Non Contactor
(Protect Glass without a lens)
- Protect Glass (no scale)
- Lens cleaning cloth
- User Manual

IDS-1000 Plus

- IDS-1000 Dermatoscope
- Non Contactor
- Scale Protect Glass
- Sleeve with a lanyard
- Lens cleaning cloth
- User Manual

	IDS-1000	IDS-1000 Plus
Polarization	Cross & Parallel	
Optics	Ø18 10X	
Field of view	18mm	
Depth of field	1.5	
Resolution	95LP	
Lighting	8 white leds	
Focus mode	fixed	
Battery	Disposable Lithium 6V 2CR5	
Smartphone & Camera Compatibility	X	O
Reticle	optional	Included
Silicone Sleeve with lanyard	optional	Included

ILLUCO

102-304 SK Ventium, #166 Gosan-ro, Gunpo-si, 15850 Gyeonggi-do, South Korea
 TEL +82 70 4922 7441 Fax +82 31 429 8826 info@illuco.co.kr www.illuco.co.kr