

MCD-A1 Dental Chair



Feature:

1. All controlled by the electric valve
2. DC Motor.(2000N)
3. Automatic thermostatic water supply system
4. Handpiece tubing with standard fittings (3sets),
5. Three way syringe (one for hot, one for cold)(2sets),
6. Water suction and saliva ejector (with switch)(1set),
7. Powerful suction apparatus (1set each)
8. Easy cleaning Integral and turn able toughened glass spittoon
9. Strong and weak dual-purpose operation light (1set),
10. X-film viewer (1 set)
11. Spirit lock tight Equilibrium device
12. Suitable for adult and children headrest.
13. Doctor stool (1set),
14. Foot pedal (1set)

Specifications:

Voltage:	110V-240V,50Hz/60Hz
Air pressure:	0.55MPa
Water pressure:	0.2MPa-0.4MPa
Highest position:	860mm
Lowest position:	440mm
Highest smallest:	Angle of dental chair (5)
Lowest smallest:	Angle of dental chair (85)
Width:	52cm
Length:	185cm max
Net weight:	200kg
Gross weight:	250kg
Packing size:	138*100*120cm

Options Accessories

- 1.High speed air turbine Handpiece (2sets)
2. Low speed air motor Handpiece. (1 set)
- 3.Fiber optic system for high speed Handpiece (1set)
4. Fiber optic Handpiece.(2sets)
- 5.Bui It-in LED Curing Light (1set)
- 6.Built-in Ultrasonic scaler (1 set)
- 7.Endoscope equipment (1 set)
- 8.Oil-less medical air compressor (1 set)
- 9.Vacuum pump for suction (1 set)
- 10.Massage function (1 set)