

DRGEM
DRGEM Corporation

ACQUIDR

DR Retrofit Solution



ACQUIDR is the digital imaging system composed of a Flat Panel Detector(FPD) and an imaging workstation with software.


The digital FPD and full-feature imaging software with excellent digital image processing, will meet all your needs in the diagnostic digital radiographic field.

Advanced Digital Image Acquisition System



- ◆ DR upgrade solution by retrofit
- ◆ Portable/fixed & wired/wireless FPD DR solution for maximum flexibility in virtually all general radiographic specialties
- ◆ The FPD fits into most existing analogue x-ray systems
- ◆ Faster workflow with the DR upgrade
- ◆ Easy interface with various x-ray generators

DR Upgrade with DRGEM Analog X-ray System



DR Interface



FPD



RadMax



Features

- Adaptable with all DRGEM analogue X-ray systems
- AED(Auto Exposure Detection) function
- Simple and convenient installation
- Intuitive graphic user interface
- Over 1300 pre-loaded exam profiles with DRGEM generator
- Configurable automatic studies
- PMS / EMR / HIS / RIS worklist support
- Multiple printing formats
- DICOM 3.0 compliant
- Brightness and contrast adjustment
- Rotation
- Horizontal and vertical flip
- Zoom and pan
- Left / right markers
- Optimized workflow
- Automatic shutters
- DAP interface support
- Exposure index (EI)

Robust and Durable Detector

Detector water & dust proof



Water proof Dust proof

High Image Quality



Product Components

Wireless FPD



Imaging Software & Workstation



*Laptop (optional)

DRGEM Customer Service Center

24/7 Support

DRGEM's Customer Service team is available 24/7, 365 days a year* via e-mail or Skype to support clients whenever they need.

*May vary by country.

Remote Diagnosis

Remote diagnosis can help minimize downtime by quickly finding the cause of errors and by detecting potential issues before problems occur.

