

REUSABLE VIDEO LARYNGOSCOPE

Stainless Steel Blades



MILLERO



MAC1



MAC2



MAC3



MAC4



D-BLADE

Blade parameter size

Blade	Mil0	MAC1	MAC2	MAC3	MAC4	D blade
Application	Neonate	Pediatric	Small Adult	Adult	Large Adult	Difficult airway
Intubation Inner Diameter	2.5-4.0mm	4.0-6.0mm	5.0-6.0mm	6.0-9.0mm	7.5-10.0mm	7.5-10.0mm
Length (cm)	8.2	10.45	12.11	14.25	16.08	13.27
Sterilize Method	2.0%Glutaral fully soak for 10 minutes Or 0.5%-0.6% Cidex OPA for more than 5 minutes					

Technical Specification

Parts	Description	Technical Index
Machine	Display	3inch LCD
	Resolution Ratio	640*480 (RGB)
	Viewing angle	≥60°
	Power	<2W
	Display front and back rotation Angle	0°-130°
	Weight	225g
Camera	Illuminance	≥150 LUX
	Resolution Ratio	1280*720 px
Image/Video function	Image/ Video function	Yes
	Output	USB output. Easy to storage and establish image/video
	Memory type	Micro SD card 8GB
Battery	Battery Type	Rechargeable lithium battery
	Capacity	1350mAh
	Cycle life of battery	>500 times
	Battery working time	>240min
	Charging time	<2(hrs)
	Charging Port	Micro USB
Power Adapter	Input	100-250V, 50 Hz.
	Output	5V1A
Working Condition	Temperature	-5 C ~+50 C
	Humidity	10%-90%
	Air Pressure	860-1060hpa
Transport / Storage Condition	Temperature	-10 C ~50 C
	Humidity	≤93%
	Air Pressure	500-1060hpa

VIDEO Laryngoscope make it better prepared for unexpected tracheal intubations.

- Confidence — MAC-style video laryngoscope blades for a wide range of clinical preferences
- Simplicity — One-button operation, Automatic power off
- Convenience — Anti-Fog Camera, 3.0inch screen reduces blind spot ,250 minutes battery with screen low battery alarming notice

Compared to traditional MAC-style laryngoscopy Blades, Video Laryngoscope articulated tip have a bigger angle that allows to raise the epiglottis equal to or greater than 60°.