


COMBISET FOR VESICAL CATHETERISATION



- Customisable pack contents and set-up
- Also available in rigid blister packs



COMBISET FOR VESICAL CATHETERISATION

COMBISET FOR VESICAL CATHETERISM

The practicality of «COMBISET» facilitates the vesical catheterization technique since it contains everything necessary for the execution of this procedure (the urological catheter is not supplied in the kit).

The sequence in which the components are found in the kit allows aseptic execution, minimising the risk of infection so easily induced in the catheterisation procedure.

The single-dose disinfectant is proprietary agent of which Bioservice is an authorised distributor and it maintains its bactericidal efficacy after sterilization.

Some versions also have a closed-circuit drainage bag which can be connected to the catheter immediately after catheterization (1).

All models can also be supplied in rigid packaging which can be used as container for the components.

The device is available in different standard configurations or is customisable according to the customer's requirements.

The device is supplied in EtO sterilized packaging.

The contents of the set do not undergo any undue alteration after sterilization.



STANDARD CONFIGURATION

CODE	MODEL	COMPOSITION	PACKAGING
H92300E	COMBISET 1	1 non-woven fabric sheet, 1 fenestrated non-woven fabric sheet with cut, 2 pairs of vinyl gloves, 2 eight-ply gauzes, 3 gauze tampons, 1 single-dose disinfectant, 1 syringe of lubricating gel, 1 syringe of distilled water, 1 forceps, 1 waste bag, 1 tray	20 pcs
H92320E	COMBISET 2	contents of "COMBISET 1" + bag model EXEL 2000-110 in transparent bag	15 pcs
H92350E	COMBISET 3	contents of "COMBISET 1" + bag model EXEL 2500-110 in transparent bag	15 pcs
H92311E	COMBISET 1 RIGID	1 non-woven fabric sheet, 1 fenestrated non-woven fabric sheet with cut, 2 pairs of vinyl gloves, 2 eight-ply gauzes, 3 gauze tampons, 1 single-dose disinfectant, 1 syringe of lubricating gel, 1 syringe of distilled water, 1 forceps, 1 waste bag	20 pcs

The individual packaging of the device consists of a rigid tray which can be used as container for the components.

CATHETERIZATION ACCESSORIES

Graduated syringes containing lubrication gel (biocompatible), water, or water and glycerine can be purchased separately as catheterization accessories.

The syringes are supplied in individual packaging sterilized by EtO.

SYRINGES

CODE	CAPACITY	CONTENTS	PACKAGING
H92301E	Syringe of 10 cc	LUBRICATING GEL	50 pcs
H92302E	Syringe of 10 cc	DISTILLED WATER	50 pcs
H92303E	Syringe of 10 cc	DISTILLED WATER WITH 10% G LYCERINE	50 pcs
H92309E	Syringe of 5 cc	LUBRICATING GEL	50 pcs
H92310E	Syringe of 5 cc	DISTILLED WATER	50 pcs
H92307E	Syringe of 5 cc	DISTILLED WATER WITH 10% G LYCERINE	50 pcs