



# BESDATA™ Dolphin

## Video Otoscopy

### Range of application:

For ear canal examination, in-hospital Ward round, out-of-hospital examination. general investigation of disease, in-hospital and out-patient consultation and surgery.

### Indications:

External auditory canal, eardrum examination and auxiliary diagnosis. Auxiliary completion of out-patient consultation, surgery.

### FEATURES:

- Easy to carry, moveable diagnosis and treatment
- Wide-angle lens, wide field of view
- Photo and video recording, easy to review
- Screen rotation for easy observation
- Replace the scope tube and quickly inserting and removing
- Special disposable devices to avoid infection
- Docking charge for easy access



## Technical Index

Display	3inch LCD
Resolution	640*480(RGB)
Power	<2W
Weight	280g
Illuminance	Natural white(6 LEDs)
Resolution	1280*720px
Image/ Video function	Yes
Output	USB output,Easy to storage and establish image/video
Memory type	Micro SD card 8GB
Battery Type	Rechargeable lithium battery
Capacity	1350mAh
Cycle life of battery	>300 times
Battery working time	>240min
Charging time	<2(hrs)

Product	Scope Model	Length of Scope(mm)	Working Length(mm)	Max OD(mm)
Video Otoscopy	VOT-01	70	30	4
Video Otoscopy	VOT-02	140	90	4.2
Video Rhinoscope	VRN-01	224	175	4.2
Video Laryngoscope	VRN-02	234	185	7