



*Creating a Better Care*

Silicone Self-adhering  
Male External Catheter  
100% Silicone

**ConSkin<sup>®</sup>**

## Greatcare offer to patients

### Security

Kink-resistant silicone funnel helps to ensure continuous urine flow, and secured connect to all major types of collector bags

### Ease of use

One piece system easy to apply and to remove, secured with an acrylic adhesive to the skin friendly.


### A perfect fit

Six different size and three styles.

Standard version for most patients and the short one for retracted penis or pediatric use.




Made of 100% silicone, the material with the highest water vapour permeability compared to latex and elasmomer



## Biocompatibility: 100% silicone

- Hypoallergenic.
- Breath ability: excellent gas exchange, less risk of maceration and irritation.
- Transparent for easy visual inspection of skin



## Measuring Devices


Measuring devices offer sizing guide to help determine the size that will ensure safe and comfortable wear.

### Standardization

CE-marking	CE
Quality system certified to	ISO 13485
Shelf life	3years

## Protective Paper Towel

Tear a small hole in the centre of a paper towel and then slip this over the penis to the base, so that hair will be not stick.



## Connect to Freedom

Made from ultrathin silicone and strong adhesive for good conformity to the penis in 3 Styles to provide greater choice and fit for a wider population.

- Style 1 (Standard) is a standard length sheath with the adhesive placed in a central band to fit to the shaft of the penis.
- Style 2 (Short Type) is a shorter sheath with the adhesive placed forward towards the tip of the sheath to suit a shorter length penis.
- Style 3 (Wide Type) is a standard sheath with a wide band of adhesive placed forward towards the tip of the penis and extending the length of the shaft.

## Adapted to all Morphologies

### ordering information

### Box of 30 Units

Style 1  
(Standard)



Size:	18mm	24mm	28mm	31mm	35mm	40mm
Ref. No.: Without Tray	GCU201461	GCU201411	GCU201412	GCU201413	GCU201414	GCU201415
Ref. No.: With Tray		GCU201416	GCU201417	GCU201418	GCU201419	GCU201420

Style 2  
(Short Type)



Size:	18mm	24mm	28mm	31mm	35mm	40mm
Ref. No.: With Tray	GCU201459	GCU201426	GCU201427	GCU201428	GCU201429	GCU201430

Style 3  
(Wide Type)



Size:	18mm	24mm	28mm	31mm	35mm	40mm
Ref. No.: With Tray	GCU201460	GCU201436	GCU201437	GCU201438	GCU201439	GCU201440





With tray




Without tray



1 pc / poly bag + paper blister



White soft blister



Transparent soft blister

### Remark:

- Standard with or without tray is available
- (Short Type) and (Wide Type) must with tray
- 1 pc of measuring guide, and 30 pcs of hair protectors are available
- OEM is available according to the customers' requirements.

