


Fiducial marker for radiation therapy



- Used to target localization for lesions or organs which move in respect to external anatomy
- Pure gold cylinders version specially knurled to inhibit migration
- Pyrolytic carbon coated zirconium oxide or tantalum barbell version to reduce the artifacts
- Provide visualization under fluoroscopy, x-ray, CT, MRI, ultrasound and other imaging modalities
- Preloaded in a pre-waxed tipped needle or in a wax-free special version
- Reduces the artifact in the MV setting (BiomarC EX)
- Ideal for proton or photon therapy (BiomarC Pro)


BiomarC GOLD




BiomarC KV ENH



BiomarC KV




BiomarC EX



BiomarC PRO



BiomarC SECURE



## BIOMARC IGRT

CODE	SIZE	DESCRIPTION	Qty/Box
040116	17Gx20cm	1,2x3mm Gold Marker	3
040118	17Gx20cm	1,2x4,2mm Gold Marker-Ergonomic handle	10
040141	18Gx20cm	0,9x3mm Gold Marker	3
040125	17Gx30cm	1,2x3mm Gold Marker	3
040117	18Gx20cm	1x5mm Carbon KV ENH Marker	3
040113	18Gx20cm	1x3mm Carbon KV Marker	3
040146	18Gx20cm	1x5mm Carbon EX Marker	3
040148	18Gx20cm	1x10mm Carbon PRO Marker	1
040150	Sutured	1,7x4,6mm Carbon Secure Marker	10

On request Vigeo is available to provide needles with different diameter and length