

# COMPACT URO

SOFT TISSUE BIOPSY

Automatic soft tissue biopsy device

CE  
0476

In compliance with MDR  
( Regulation EU 2017/745 )



- One-handed design simplifies loading and operation whilst allowing for convenient manipulation of the ultrasound transducer
- Thin-wall cannula for larger core
- Colour coded for quick gauge size identification
- Depth markings etched on cannula for precise positioning
- Loaded status indicator
- Plastic safety protection to avoid accidental firing



Flat design for transrectal ultrasound (TRUS)-guided needle biopsy



Colour code

## COMPACT URO Qty/Box 10pcs

CODE	SIZE
● VCTA1620	16Gx20cm
● VCTA1625	16Gx25cm
● VCTA1820	18Gx20cm
● VCTA1825	18Gx25cm
● VCTA1830	18Gx30cm

On request Vigeo is available to provide needles with different diameter and length