

Chiba needle for aspiration biopsy



- Puncture cannula developed especially for ultrasonic diagnostics
- Echogenic needle tip assists precise needle tip positioning
- Safe for patient and user, thanks to the adjustable sliding depth stopper
- Graduated puncture cannula to assist in precise depth placement of the needle
- Transparent luer cone for visual inspection of aspirate
- Colour coded hub for easy identification of needle gauge size



Reusable syringe gun (code: SG01)

CYTOCAN Qty/Box 20pcs

CODE	SIZE
● VVB1810	18Gx10cm
● VVB1815	18Gx15cm
● VVB1820	18Gx20cm
● VVB2010	20Gx10cm
● VVB2015	20Gx15cm
● VVB2020	20Gx20cm
● VVB2105	21Gx5cm
● VVB2110	21Gx10cm
● VVB2115	21Gx15cm
● VVB2120	21Gx20cm
● VVB2210	22Gx10cm
● VVB2215	22Gx15cm
● VVB2220	22Gx20cm
● VVB2310	23Gx10cm
● VVB2315	23Gx15cm
● VVB2505	25Gx5cm

On request Vigeo is available to provide needles with different diameter and length