

One-handed Menghini Aspiration Device



- One-handed safe and fast operation
- Pyramidal tip allows for easy insertion
- Menghini design cannula for cylindrical and large sampling
- Intended for use in liver, lymph nodes and thyroid
- Colour coded for quick gauge size identification
- Centimeter depth markings assist in needle placement



One-handed design simplifies loading and operation whilst allowing for convenient manipulation of the ultrasound transducer

HEPASHOT Qty/Box 10pcs

CODE	SIZE
● VVC01610	16Gx10cm
● VVC01615	16Gx15cm
● VVC01620	16Gx20cm
● VVC01805	18Gx5cm
● VVC01810	18Gx10cm
● VVC01815	18Gx15cm
● VVC01820	18Gx20cm
● VVC02005	20Gx5cm
● VVC02010	20Gx10cm
● VVC02015	20Gx15cm
● VVC02020	20Gx20cm
● VVC02105	21Gx5cm
● VVC02110	21Gx10cm
● VVC02115	21Gx15cm

On request Vigeo is available to provide needles with different diameter and length