

INTRO/INTRO PLUS SOFT TISSUE BIOPSY

CE
0476

In compliance with MDR
(Regulation EU 2017/745)

Standard coaxial introducer needle used for single access in multiple biopsies



- Echogenic treatment for ultrasound guided procedures
- Centimeter marking and optional sliding stopper
- Sterilized by ETO, shelf life 5 years
- INTRO (code VVT) with blunt cannula and trocar tip
- INTRO PLUS (code VVTC) with sharp cutting cannula and trocar tip



Depth marker with code colour

INTRO/INTRO PLUS Qty/Box 20pcs

CODE	SIZE
● VVT1105 VVTC1105	11Gx5cm
● VVT1305 VVTC1305	13Gx5cm
● VVT1310 VVTC1310	13Gx10cm
● VVT1315 VVTC1315	13Gx15cm
● VVT1505 VVTC1505	15Gx5cm
● VVT1510 VVTC1510	15Gx10cm
● VVT1515 VVTC1515	15Gx15cm
● VVT1705 VVTC1705	17Gx5cm
● VVT1710 VVTC1710	17Gx10cm
● VVT1715 VVTC1715	17Gx15cm

On request Vigeo is available to provide needles with different diameter and length