


Chiba needles for prenatal diagnosis



- Individually packaged in sterile wrappers
- Intended for use in amniocentesis or chorionic villus sampling
- Safe for patient and user, thanks to the adjustable sliding depth stopper
- Graduated puncture cannula to assist in precise depth placement of the needle
- Transparent luer cone for visual inspection of aspirate
- Colour coded hub for easy identification of needle gauge size
- Sharp, non-coring needle tip to minimize patient trauma
- Self-locating stylets to prevent blockage

## AMNIO



## CHORIAL KIT



## PRENATAL KIT



### AMNIO Qty/Box 20pcs

CODE	SIZE
● VVBP1915	19Gx15cm
● VVBP1920	19Gx20cm
● VVBP2015	20Gx15cm
● VVBP2020	20Gx20cm
● VVBP2115	21Gx15cm
● VVBP2120	21Gx20cm
● VVBP2315	23Gx15cm

### CHORIAL KIT Qty/Box 20pcs

● VVBU1815	18Gx15cm
------------	----------

### PRENATAL KIT Qty/Box 20pcs

● VVBL2016	20Gx16cm
● VVBL2020	20Gx20cm
● VVBL2116	21Gx16cm
● VVBL2120	21Gx20cm

On request Vigeo is available to provide needles with different diameter and length