

# RELOCK PREMIUM SOFT WOMEN'S HEALTH

CE  
0476




In compliance with MDR  
( Regulation EU 2017/745 )

PATENT PENDING

Preoperative repositionable breast localization device



- Preoperative marking of local lesions
- One-handed design for a convenient manipulation of the ultrasound transducer
- Retractable hook with unique flexibility of its polyamide guiding thread
- Indicator on wing shows direction of hook
- Echogenic needle tip for ultrasound use
- Hook release at the needle tip for accurate targeting
- Centimeter depth markings assist in needle placement
- Patented stopper for partial hook release prior to confirming position
- 9mm loop hook
- Single packaged stiffener cannula used in the operating room



Mammography

## RELOCK PREMIUM SOFT Qty/Box 10pcs

CODE	SIZE
● VVRPS1905	19Gx5cm
● VVRPS1907	19Gx7cm
● VVRPS1910	19Gx10cm
● VVRPS2005	20Gx5cm
● VVRPS2007	20Gx7cm
● VVRPS2010	20Gx10cm

On request Vigeo is available to provide needles with different diameter and length

1. Half exposure - 2. Full exposure - 3. Stiffener cannula