



## Preoperative repositionable breast localization device



- Preoperative marking of local lesions
- Retractable hook with unique flexibility
- Indicator shows direction of hook
- Echogenic needle tip for ultrasound use
- Hook release at the needle tip for accurate targeting
- Centimeter depth markings assist in needle placement
- Wire fixation clamp to avoid migration
- 9mm loop hook



Mammography



Wire fixation clamp

### RELOCK Qty/Box 10pcs

CODE	SIZE
● VVR1905	19Gx5cm
● VVR1907	19Gx7cm
● VVR1910	19Gx10cm
● VVR2007	20Gx7cm
● VVR2010	20Gx10cm
● VVR2015	20Gx15cm

On request Vigeo is available to provide needles with different diameter and length