

One-handed preoperative breast localization device



- Preoperative marking of local lesions
- One-handed design for a convenient manipulation of the ultrasound transducer
- Puncture cannula with ultrasonic tip with guidance wire in place
- The Z-shaped hook provides safer delivery in critical areas (available in 18G version)
- The X-shaped hook provides excellent fixation and is intended to limit migration (available in 20 and 21G version)



Single Hook



Z-shaped hook



X-shaped hook



One-handed design

SURGIMARC SLIM <small>Qty/Box 10pcs</small>		
CODE		SIZE
● VVNS1807	VVNSX1807 (Z-shaped hook)	18Gx7cm
● VVNS1810	VVNSX1810 (Z-shaped hook)	18Gx10cm
● VVNS2007	VVNSX2007 (X-shaped hook)	20Gx7cm
● VVNS2010	VVNSX2010 (X-shaped hook)	20Gx10cm
● VVNS2015	VVNSX2015 (X-shaped hook)	20Gx15cm
● VVNS2107	VVNSX2107 (X-shaped hook)	21Gx7cm
● VVNS2110	VVNSX2110 (X-shaped hook)	21Gx10cm

On request Vigeo is available to provide needles with different diameter and length