



- Safety button to prevent unwanted shots
- Two triggers
- Two possible puncture depths of 16 or 22 mm
- Compatible exclusively with Vigeo Autocut biopsy needle
- Loaded status indicator



Loaded status indicator



Coaxial kit (VGC)

| VTEK | | |
|------|----------|--------|
| CODE | THROW | SAFETY |
| VTEK | 16-22 mm | YES |

| AUTOCUT Qty/Box 10pcs | | |
|-----------------------|----------|-------------|
| CODE | SIZE | COAXIAL KIT |
| ● VVG1210 | 12Gx10cm | VVGC1210 |
| ● VVG1410 | 14Gx10cm | VVGC1410 |
| ● VVG1413 | 14Gx13cm | VVGC1413 |
| ● VVG1415 | 14Gx15cm | VVGC1415 |
| ● VVG1610 | 16Gx10cm | VVGC1610 |
| ● VVG1615 | 16Gx15cm | VVGC1615 |
| ● VVG1620 | 16Gx20cm | VVGC1620 |
| ● VVG1810 | 18Gx10cm | VVGC1810 |
| ● VVG1815 | 18Gx15cm | VVGC1815 |
| ● VVG1820 | 18Gx20cm | VVGC1820 |
| ● VVG1825 | 18Gx25cm | VVGC1825 |
| ● VVG2010 | 20Gx10cm | VVGC2010 |
| ● VVG2015 | 20Gx15cm | VVGC2015 |

On request Vigeo is available to provide needles with different diameter and length