

PHYSIOLASER 12.0 PRO

LASER THERAPY


THE "PRO" VERSION OF PHYSIOLASER 12.0 IS DESIGNED FOR THE PROFESSIONAL WHO USES THE DEVICE FREQUENTLY AND PREFERS TO KEEP IT ALWAYS NEARBY ON A CART. THE 12W OF MAXIMUM POWER ALLOW RAPID AND EFFECTIVE TREATMENTS EVEN ON THE LARGEST BODY AREAS, USING THE 22 MM OPTIC INCLUDED. THE WIDE LIBRARY OF PROGRAMS ALSO GUARANTEES THE PROFESSIONAL TO TREAT ANY MUSCULOSKELETAL DISORDER AND SOME PROBLEMS RELATING TO WOUNDS AND ULCERS OF THE SKIN.

TECHNICAL SPECIFICATIONS

Display
2.6", backlit
Laser source
Class 4 GaAlAs diode
Maximum power
12W ±20%
Compatible with a wavelength of
808nm 980nm and 1064nm (optional)
Issue
Pulsed and continuous
Target area dimensions
11 mm optic 22 mm optic
Frequency
Up to 10,000Hz
Power supply
Mains and rechargeable batteries
Size and weight
220x170x60 mm - 1.12 Kg

STANDARD EQUIPMENT

- 1 Suitcase
- 1 PHYSIOLASER 12.0 PRO device
- 1 808nm laser handpiece with 11mm and 22mm optics
- 2 Protection glasses (patient and operator)
- 1 Emergency button
- 1 Power supply unit
- 1 Operating manual



PURCHASE OPTIONS

- Code G5995 HP 980nm handpiece
- Code G5996 1064nm handpiece, maximum power 3W



PHYSIOLASER 12.0 PRO

LASER THERAPY

MEDICAL PROGRAM LIST

808nm WAVELENGTH PROGRAMS

MUSCULOSKELETAL

GENERAL

Arthrosis
Rheumatoid arthritis
Muscle injuries
Muscle contracture
Decontracting
Fibromyalgia
Radiculopathies
Tendinopathies
Neuralgia
Postherpetic neuralgia
Sudeck disease
Bursitis
Edema

ARM - HAND

Arthrosis
Rheumatoid arthritis
Raynaud syndrome
Finger pain
Epicondylitis
Carpal tunnel
Wrist fractures
Wrist pain
Tendinitis
Edema
Radius fracture

KNEE

Arthrosis
Rheumatoid arthritis
Acute pain
Chronic pain
Meniscus pain
Bursitis
Patello-femoral
Edema
Ligament injury
Gonalgia
Sprain
Chondromalacia
Chondropathy
Patello-femoral tendinitis

FOOT - ANKLE

Arthrosis
Rheumatoid arthritis
Metatarsalgia
Ankle pain
Plantar fasciitis
Achilles tendon
Sprain
Edema
Heel pain

BACK

Vertebral arthrosis
Acute pain
Lumbago
Myofascial pain
Radiculopathies
Rheumatoid arthritis

SHOULDER

Myofascial pain
Arthrosis
Adhesive capsulitis
Bursitis
Subacromial impingement
Tendinitis
Luxations
Tendon injuries (cuff)
Rheumatoid arthritis
Edema

HEAD - NECK

Arthrosis
Radiculopathies
Myofascial pain
Chronic cervical
Temporo-mandibular pain
Rheumatoid arthritis
Trigeminal neuralgia
Cervical herniation

DERMATOLOGY

Ulcers
Cicatrizazion ulcers
Wounds

980nm WAVELENGTH PROGRAMS

Antalgic
Cervical Pain
Lumbar pain
Bursitis
Baker's cyst
Cruralgia
De Quervain syndrome
Sprains
Contusions
Acute pain
Chronic pain
Edemas
Hematomas
Extended hematomas
Epicondylitis
Epitrochleitis
Discal hernia
Fibromyalgia
Inflammations
Ligament injury
Muscle injuries
Meniscus pain
Synovitis
Muscle strains
Tendinopathy
Carpal Tunnel