

PHYSIOLASER 2.0

LASER THERAPY

LIGHT AND PORTABLE DEVICE, EQUIPPED WITH ALL THE MAIN PROGRAMS TO TREAT CHRONIC AND ACUTE PATHOLOGIES OF THE MUSCULOSKELETAL SYSTEM.

TECHNICAL SPECIFICATIONS

Display

2,6", backlit

Laser source

Class 4 GaAlAs diode

Maximum power

2W \pm 20%

Compatible with a wavelength of

808nm

Issue

Pulsed and continuous

Target area dimensions

11 mm optic

Frequency

Up to 10,000Hz

Power supply



Mains and rechargeable batteries

Size and weight

100x160x35 mm - 450 g

STANDARD EQUIPMENT

- 1 Suitcase
- 1 PHYSIOLASER 2.0 device
- 1 808nm laser handpiece with 11mm optic
- 2 Protection glasses (patient and operator)
- 1 Emergency button
- 1 Power supply unit
- 1 Operating manual



CE
0476



FAVORITES



PROGRAMMABLE



BATTERY PACK



MAINS POWER SUPPLY

GLOBUS
ITALIAN EXCELLENCE

PHYSIOLASER 2.0

LASER THERAPY

MEDICAL PROGRAM LIST

MUSCULOSKELETAL

GENERAL

Arthrosis

Rheumatoid arthritis

Muscle contracture

Fibromyalgi

Tendinopathies

Bursitis

Edema

ARM - HAND

Arthrosis

Rheumatoid arthritis

Epicondylitis

Carpal tunnel

Tendinitis

KNEE

Arthrosis

Rheumatoid arthritis

Meniscopathy

Patello-femor

Ligament injury

Gonalgia

Chondropathy

FOOT - ANKLE

Arthrosis

Rheumatoid arthrit

Ankle pain

Plantar fasciitis

Achilles tendinitis

Sprain

BACK

Vertebral arthrosis

Acute pain

Lumbago

Myofascial pain

Rheumatoid arthritis

SHOULDER

Myofascial pain

Arthrosis

Adhesive capsulitis

Subacromial conflict

Tendinitis

Rheumatoid arthriti

HEAD - NECK

Arthrosis

Arthrosis

Rheumatoid arthritis

DERMATOLOGY

Wounds