

# PHYSIOLASER 500

## LASER THERAPY

THE PRODUCT IS INDICATED TO ACCELERATE THE RECOVERY OF INFLAMED TISSUES, REGENERATE CELLS AND IMPROVE MICROCIRCULATION. PHYSIOLASER 500 PROMOTES THE NATURAL HEALING PROCESSES.

### TECHNICAL SPECIFICATIONS

#### Display

2,6", backlit

#### Laser source

Class 4 GaAlAs diode

#### Maximum power

0,5W ±20%

#### Compatible with a wavelength of

808nm

#### Issue

Pulsed and continuous

#### Target area dimensions

11 mm optic

#### Frequency

Up to 10,000Hz

#### Power supply


Mains and rechargeable batteries

#### Size and weight

100x160x35 mm - 450 g

### STANDARD EQUIPMENT

- 1 Suitcase
- 1 PHYSIOLASER 500 device
- 1 808nm laser handpiece with 11mm optic
- 2 Protection glasses (patient and operator)
- 1 Emergency button
- 1 Power supply unit
- 1 Operating manual



CE  
0476



FAVORITES



PROGRAMMABLE



BATTERY PACK



MAINS POWER SUPPLY

**GLOBUS**  
ITALIAN EXCELLENCE

# PHYSIOLASER 500

## LASER THERAPY

### PROGRAM LIST

#### HEALTH - REHAB

Arthrosis

---

Ankle sprain

---

Knee pain

---

Edema Ankle sprain

---

Shoulder articular pain

---

Epicondylitis

---

Chronic pain Back pain

---

Neck pain

---

Tendinopathies

#### AESTHETIC - FITNESS

##### FACE AESTHETIC

Wrinkles

---

Acne

---

Skin rejuvenation

---

Skin - face lift

---

Expression wrinkles

---

Rosacea

##### BODY AESTHETIC

Breast stretch marks

---

Thighs, hips and gluteus edematous cellulite

---

Non-localized body fat and cellulite

---

Stretch marks

---

Scars