

PHYSIOLASER 6.0 PRO

LASER THERAPY


PHYSIOLASER 6.0 PRO IS A DEVICE OF SMALL WEIGHT AND SIZE, SUITABLE BOTH TO BE PLACED ON A STUDIO CART AND TO BE TRANSPORTED, IF NECESSARY, TO OTHER LOCATIONS. WITH ITS MAXIMUM POWER OF 6W, IT IS IDEAL FOR TREATING OSTEOARTICULAR AND MUSCULAR PATHOLOGIES, QUICKLY REDUCING PAIN AND EDEMA AND ACCELERATING THE RECOVERY OF INJURED AND INFLAMED TISSUES.

TECHNICAL SPECIFICATIONS

Display
2.6", backlit
Laser source
Class 4 GaAlAs diode
Maximum power
6W \pm 20%
Compatible with a wavelength of
808nm 980nm and 1064nm (optional)
Issue
Pulsed and continuous
Target area dimensions
11 mm optic 22 mm optic (optional)
Frequency
Up to 10,000Hz
Power supply
Mains and rechargeable batteries
Size and weight
220x170x60 mm - 1.12 Kg

STANDARD EQUIPMENT



- 1 Suitcase
- 1 PHYSIOLASER 6.0 PRO device
- 1 808nm laser handpiece with 11mm optic
- 2 Protection glasses (patient and operator)
- 1 Emergency button
- 1 Power supply unit
- 1 Operating manual



PURCHASE OPTIONS

- Cod. G5994 MP 980nm handpiece
- Cod. G5996 1064nm handpiece, maximum power 3W
- Cod. G6092 22 mm optic



-  FAVORITES
-  PROGRAMMABLE
-  BATTERY PACK
-  MAINS POWER SUPPLY

PHYSIOLASER 6.0 PRO

LASER THERAPY

MEDICAL PROGRAM LIST

808nm WAVELENGTH PROGRAMS

MUSCULOSKELETAL

GENERAL

Arthrosis

Rheumatoid arthritis

Muscle injuries

Muscle contracture

Radiculopathies

Tendinopathies

Neuralgia

Bursitis

Edema

ARM - HAND

Arthrosis

Rheumatoid arthritis

Epicondylitis

Carpal tunnel

Wrist pain

Tendinitis

Edema

KNEE

Arthrosis

Rheumatoid arthritis

Acute pain

Chronic pain

Meniscus pain

Bursitis

Patello-femoral

Edema

Ligament injury

Gonalgia

Sprain

Chondropathy

Patello-femoral tendinitis

FOOT - ANKLE

Arthrosis

Rheumatoid arthritis

Ankle pain

Plantar fasciitis

Achilles tendon

Sprain

Edema

Heel pain

BACK

Vertebral arthrosis

Acute pain

Lumbago

Myofascial pain

Radiculopathies

Rheumatoid arthritis

SHOULDER

Myofascial pain

Arthrosis

Adhesive capsulitis

Bursitis

Subacromial impingement

Tendinitis

Tendon injuries (cuff)

Rheumatoid arthritis

HEAD - NECK

Arthrosis

Radiculopathies

Chronic cervical

Temporo-mandibular pain

Rheumatoid arthritis

Trigeminal neuralgia

DERMATOLOGY

Ulcers

Wounds

980nm WAVELENGTH PROGRAMS

Antalgic

Cervical Pain

Lumbar pain

Bursitis

Baker's cyst

Cruralgia

De Quervain syndrome

Sprains

Contusions

Acute pain

Chronic pain

Edemas

Hematomas

Extended hematomas

Epicondylitis

Epitrochleitis

Discal hernia

Fibromyalgia

Inflammations

Ligament injury

Muscle injuries

Meniscus pain

Synovitis

Muscle strains

Tendinopathy

Carpal Tunnel