

PODCARE 2.0

LASER THERAPY


PODCARE 2.0 IS A PORTABLE DEVICE WITH PROTOCOLS TO TREAT DIFFERENT MUSCULOSKELETAL PATHOLOGIES AND WITH SPECIAL PODOLOGICAL APPLICATIONS SUCH AS ONICOMYCOSIS, ONICOCRIPTOSIS AND OTHER TYPICAL INFLAMMATORY PATHOLOGIES OF THE FOOT.

TECHNICAL SPECIFICATIONS

Display
2.6", backlit
Laser source
Class 4 GaAlAs diode
Maximum power
2W ±20%
Compatible with a wavelength of
808nm
Issue
Pulsed and continuous
Target area dimensions
11 mm optic 3 mm optic (optional)
Frequency
Up to 10,000Hz
Power supply
Mains and rechargeable batteries
Size and weight
100x160x35 mm - 450 g

STANDARD EQUIPMENT

- 1 Suitcase
- 1 PODCARE 2.0 device
- 1 808nm laser handpiece with 11mm optic
- 2 Protection glasses (patient and operator)
- 1 Emergency button
- 1 Power supply unit
- 1 Operating manual



PURCHASE OPTIONS

- Cod. G5547 3 mm optic



CE
0476



FAVORITES



PROGRAMMABLE



BATTERY PACK



MAINS POWER SUPPLY

GLOBUS
ITALIAN EXCELLENCE

PODCARE 2.0

LASER THERAPY

MEDICAL PROGRAM LIST

MUSCULOSKEL

GENERAL

Arthrosis

Rheumatoid arthritis

Muscle contracture

Fibromyalgia

Radiculopathies

Tendinopathies

Neuralgia

Bursitis

Edema

ARM - HAND

Arthrosis

Rheumatoid arthritis

Epicondylitis

Carpal tunnel

Wrist pai

KNEE

Arthrosis

Rheumatoid arthritis

Meniscopathy

Ligament injury

Chondropathy

FOOT - ANKLE

Arthrosis

Rheumatoid arthrit

Onychocryptosis

Onychomycosis

Ankle pain

Plantar fasciitis

Achilles tendinitis

Sprain

Morton's neuroma

BACK

Vertebral arthrosis

Lumbago

Radiculopathies

SHOULDER

Adhesive capsulitis

Subacromial impingement

Tendon lesions (cuff)

Rheumatoid arthritis

HEAD - NECK

Arthrosis

Chronic cervical

Temporo-mandibular pain

Rheumatoid arthritis

DERMATOLOGY

Wounds