

# MAGNETOVET 4000

## MAGNETOTHERAPY

MAGNETOVET 4000 IS A HIGH INTENSITY PULSED ELECTROMAGNETIC FIELD MAGNETOTHERAPY EQUIPMENT. IT IS POSSIBLE TO USE FOUR OUTPUT CHANNELS ALLOWING TO CONNECT MANY ACCESSORIES AND TO TREAT MORE BODY PARTS AT THE SAME TIME OR EVEN TWO PETS SIMULTANEOUSLY. MAGNETOTHERAPY HAS COUNTLESS AND PROVEN BENEFICIAL EFFECTS ON BOTH BONE AND ARTICULAR TISSUE. ELECTROMAGNETIC IMPULSES PRODUCED, IN FACT, STIMULATE THE REGENERATION OF NEW CELLS AND THE REMODELING OF EXISTING STRUCTURES.


SPECIFIC PROGRAMS FOR DOGS, CATS AND HORSES

### TECHNICAL FEATURES

<b>Display</b>
Backlit, 2.6"
<b>Channels</b>
4
<b>Emitted intensity</b>
200 average Gauss per channel with flex solenoids 800 total average Gauss (for the 4 channels)
<b>Frequency</b>
Adjustable in programming from 5 to 200Hz
<b>Duty Cycle</b>
From 5% to 50% in steps of 5
<b>Power supply</b>
- Mains - Battery pack (optional)
<b>Size</b>
170x220x60 mm
<b>Weight</b>
1370 gr


### STANDARD EQUIPMENT

- 1 Carrying case
- 1 MAGNETOVET 4000 device
- 4 Flexible solenoids 300x100 mm
- 1 Power supply unit
- 1 Operating manual



### PURCHASE OPTIONS

- Cod. G6039 Memory-foam therapeutic bed, 4 coils, 700x440 mm, for cats and small dogs.. For large dogs, buy 2 single beds and put them together with the velcro.



- LAST 10
- MULTI USER
- PROGRAMMABLE
- RENTAL OPTION
- 2+2 MODE
- MAINS POWER SUPPLY

# MAGNETOVET 4000

## MAGNETOTHERAPY

Rev. En\_03\_2021

### PROGRAM LIST

#### DOG / CAT

Fractures

Pseudoarthrosis

Osteoarthrosis

Wounds / Ulcers

Hip osteoarthrosis

Knee osteoarthrosis

Elbow osteoarthrosis

Shoulder osteoarthrosis

Spondylarthrosis

Discopathy

Luxation

Arthritis

Arthrosis

Delayed Fracture Consolidation

Cicatrization

Edema

Muscle relaxation

Urinary retention

Smooth muscles

Peripheral paralyses

Inflammation pain

#### HORSE

Fractures

Wounds / Ulcers

Elbow osteoarthrosis

Shoulder osteoarthrosis

Arthritis

Arthrosis

Delayed Fracture Consolidation