

TECARVET 2000

TECAR THERAPY

TECARVET 2000 COMBINES ALL THE QUALITY OF GLOBUS TECAR THERAPY IN A MACHINE OF REDUCED DIMENSIONS AND WEIGHT, TO GUARANTEE MAXIMUM PORTABILITY, KEEPING ITS EFFECTIVENESS AND PROFESSIONALISM UNCHANGED. THE SPECIAL RESISTIVE BIPOLAR HANDPIECE ALLOWS TO SPEED UP PREPARATION FOR TREATMENT.

SPECIFIC PROGRAMS FOR DOGS AND CATS

TECHNICAL FEATURES

Display

Backlit, 2.6"

Power

50W, $\pm 10\%$

Emission Frequency

470KHZ, $\pm 10\%$

Treatment modality

Resistive bipolar

Power supply

- Mainse
- Battery pack

Size


170x220x60 mm

Weight

1100 gr

STANDARD EQUIPMENT

- 1 Professional Carrying case
- 1 TECARVET 2000 device
- 1 bipolar handpiece \varnothing 50 mm
- 1 260 ml gel bottle
- 1 Power supply unit
- 1 Operating manual



CE



-  RESISTIVE
-  PROGRAMMABLE
-  BATTERY PACK
-  MAINS POWER SUPPLY

TECARVET 2000

TECAR THERAPY

PROGRAM LIST

DOG / CAT

Osteoarthritis

Muscle contractures

endinopathies

Edema