

ELITE

4 CHANNELS ELECTROTHERAPY | 98 PROGRAMS

FOUR CHANNELS PRODUCT FOR EVERY TREATMENT NEEDS.

32 SPORT | 29 BEAUTY | 26 FITNESS | 9 TENS | 2 REHAB

TECHNICAL FEATURES

Display

Visible area size 2,6"

Channels

4 independent (8 electrodes)

Frequency

0,3-150 Hz *

Pulse amplitude

40-450 μ s *

Power

0-120 mA per channel

Power supply

Rechargeable batteries

Size

100x160x35 mm


Weight

440 gr

* According to the type of current

EQUIPMENT

- 1 Bag
- 1 Stimulator ELITE
- 4 Cables for electrode connection
- 4 Self-adhesive square electrodes
- 4 Self-adhesive rectangular electrodes
- 1 Charger
- 1 Operating manual



CE
0476



- TENS
- LAST 10
- MULTI USER
- MY TRAINER
- SYNCRO STIM
- RECHARGEABLE BATTERIES

GLOBUS
ITALIAN EXCELLENCE

PROGRAMS LIST

SPORT

- Demo
- Motor point pen
- Capillarization
- Warm-up
- Decontracting
- Post-competition recovery/training
- Maximum strength
- Endurance strength
- Explosive strenght

Reactivity

FITNESS

- Firming
- Bio-Pulse Firming
- Sculpting
- Toning
- Mass building

BEAUTY

- Drainage
- Bio-Pulse Drainage
- Lipolysis
- Bio-Pulse Lipolysis
- Toning massage
- Bio-Pulse relaxing massage
- Vascularization
- Face Capillarization
- Face lifting effect

TENS AND PAIN

- Antalgic tens
- Endorphinic tens
- Contractures
- Hematomas
- Sciatica
- Lumbago
- Cervical pain
- Shoulder pain
- Epicondylitis

REHAB

- Lower limbs atrophy recovery
- Upper limbs atrophy recovery