

# LASERVET 1.0

## LASER THERAPY

LASERVET 1.0 IS THE ONLY LASER DEVICE IN CLASS 4 OF THE GLOBUSVET CATALOGUE. THE DEVICE CAN REACH 1 WATT OF POWER AND IT INCLUDES PROGRAMS FOR DOGS, CATS AND HORSES. MOREOVER, IT IS POSSIBLE TO CREATE AND SAVE NEW PROGRAMS.

SPECIFIC PROGRAMS FOR DOGS, CATS AND HORSES

### TECHNICAL FEATURES

#### Display

Backlit, 2.6"

#### Laser source

Class 4 GaAlAs diode

#### Peak power

1W  $\pm$ 20%

#### Preset for wave length

808nm

#### Issue

Pulsed and continuous

#### Target area dimensions

11 mm  $\varnothing$ , optics

#### Frequency

Up to 10000Hz

#### Power supply

Mains and rechargeable battery

#### Size


100x160x35 mm

#### Weight

450 gr







### STANDARD EQUIPMENT

- 1 Carrying case
- 1 LASERVET 1.0 device
- 1 Handpiece White 808nm with 11 mm  $\varnothing$ , optics
- 2 Protection glasses
- 1 Power supply unit
- 1 Operating manual



CE



-  LAST 10
-  MULTI USER
-  PROGRAMMABLE
-  RENTAL OPTION
-  BATTERY PACK
-  MAINS POWER SUPPLY

# LASERVET 1.0

## LASER THERAPY

**MORE THAN 40 PROGRAMS FOR DOGS, CATS, HORSES.**  
**TREATMENTS INCLUDE PROTOCOLS FOR:**

**Muscle-skeletal injuries** such as arthrosis, tendinopathies, contusions, bursitis, sprains, muscle strains, etc.

**Pathologies of the integumentary system** such as skin or deep wounds

**Pathologies of other systems** such as otitis, sinusitis, stomatitis etc.

