

LASERVET 3.0

LASER THERAPY

THE TREATMENT WITH LASERVET 3.0 ACCELERATES THE HEALING OF THE DAMAGED TISSUE AND IMPROVES CIRCULATION. LASERVET 3.0 CAN BE USED TO FACILITATE THE CICATRIZATION OF SKIN WOUNDS AND IT IS PARTICULARLY EFFECTIVE IN THE TREATMENT OF DIFFERENT MUSCULOSKELETAL PATHOLOGIES CHARACTERIZED BY INFLAMMATION AND PAIN. LASERTHERAPY ACCELERATES THE HEALING AND THE PAIN ELIMINATION WHEN USED ALONE AND ALSO TOGETHER WITH OTHER THERAPEUTIC MODALITIES.


SPECIFIC PROGRAMS FOR DOGS, CATS, HORSES, SMALL MAMMALS, REPTILES AND BIRDS

TECHNICAL FEATURES

Display
Backlit, 2.6"
Laser source
Class 4 GaAlAs diode
Peak power
3W ±20%
Preset for wave length
808nm
Issue
Pulsed and continuous
Target area dimensions
11 mm Ø, optics
Frequency
Up to 10000Hz
Power supply
Mains and rechargeable battery
Size
100x160x35 mm
Weight
450 gr

STANDARD EQUIPMENT

- 1 Carrying case
- 1 LASERVET 3. device
- 1 Handpiece Black 808nm with 11 mm Ø, optics
- 2 Protection glasses
- 1 Power supply unit
- 1 Operating manual



- LAST 10
- MULTI USER
- PROGRAMMABLE
- RENTAL OPTION
- BATTERY PACK
- MAINS POWER SUPPLY

LASERVET 3.0

LASER THERAPY

MORE THAN 130 PROGRAMS FOR DOGS, CATS, HORSES, SMALL MAMMALS, REPTILES AND BIRDS. TREATMENTS INCLUDE PROTOCOLS FOR:

Muscle-skeletal injuries such as arthrosis, tendinopathies, contusions, bursitis, sprains, muscle strains, fractures, dysplasia, neuropathies etc.

Pathologies of the integumentary system such as skin wounds, granulomas, sacculitis, perianal fistulas, pododermatitis, abrasions, etc.

Pathologies of other systems such as otitis, sinusitis, stomatitis etc.

