

# LASERVET 8.0

## LASER THERAPY

LASERVET 8.0 IS A PORTABLE AND HANDHELD DEVICE FOR HIGH PULSED POWER LASER THERAPY THAT ALLOWS TO PERFORM VETERINARY TREATMENTS. THE ACTION OF LASER IS AIMED AT TISSUE REGENERATION, POST-TRAUMATIC AND ANTI-INFLAMMATORY TREATMENT AND REDUCTION OF BOTH CHRONIC AND ACUTE PAIN.

SPECIFIC PROGRAMS FOR DOGS, CATS, HORSES, SMALL MAMMALS, REPTILES AND BIRDS

### TECHNICAL FEATURES

#### Display

Backlit, 2.6"

#### Laser source

Class 4 GaAlAs diode

#### Peak power

8W  $\pm$ 20%

#### Preset for wave length

808nm

980nm and 1064nm (optional)

#### Issue

Pulsed and continuous

#### Target area dimensions

11 mm  $\varnothing$ , optics

22 mm  $\varnothing$  and 60 mm  $\varnothing$ , optics (optional)

#### Frequency

Up to 10000Hz

#### Power supply

Mains and rechargeable battery

#### Size


100x160x35 mm

#### Weight

450 gr

### STANDARD EQUIPMENT

- 1 Carrying case
- 1 LASERVET 8.0 device
- 1 Handpiece 808nm with 11 mm  $\varnothing$ , optics
- 2 Protection glasses
- 1 Power supply unit
- 1 Operating manual









### PURCHASE OPTIONS

- Cod. G6347 Handpiece MP 980nm
- Cod. G5996 Handpiece 1064nm
- Cod. G6089 Optics 60 mm  $\varnothing$
- Cod. G6092 Optics 22 mm  $\varnothing$

CE



-  LAST 10
-  MULTI USER
-  PROGRAMMABLE
-  RENTAL OPTION
-  BATTERY PACK
-  MAINS POWER SUPPLY

# LASERVET 8.0

## LASER THERAPY

**MORE THAN 130 PROGRAMS FOR DOGS, CATS, HORSES, SMALL MAMMALS, REPTILES AND BIRDS. TREATMENTS INCLUDE PROTOCOLS FOR:**

**Muscle-skeletal injuries** such as arthrosis, tendinopathies, contusions, bursitis, sprains, muscle strains, fractures, dysplasia, neuropathies etc.

**Pathologies of the integumentary system** such as skin wounds, granulomas, sacculitis, perianal fistulas, pododermatitis, abrasions, etc.

**Pathologies of other systems** such as otitis, sinusitis, stomatitis etc.

