

DUO TENS

ELECTROTHERAPY 2 CHANNELS | 20 PROGRAMS

THE IDEAL PRODUCT TO TREAT BACK PAIN AND CERVICAL PAIN

16 TENS AND PAIN | 3 SPORT | 1 FITNESS, BEUTY, WELLNESS

TECHNICAL FEATURES

Display

Visible area size 2,5"

Channels

2 independents (4 electrodes)

Frequency

1-150Hz

Pulse amplitude

50-400µs

Power

0-100mA per channel

Power supply

Rechargeable battery

Size


mm 80x130x25

Weight

gr 220





EQUIPMENT

- 1 Bag
- 1 DUO TENS stimulator
- 2 Cables for electrode connection
- 4 self-adhesive square electrodes
- 4 self-adhesive rectangular electrodes
- 1 Charger
- 1 Operating manual



CE
0476



-  TENS
-  PROGRAMMABLE
-  SYNCRO STIM
-  RECHARGEABLE BATTERY

GLOBUS
ITALIAN EXCELLENCE

DUO TENS

ELECTROTHERAPY 2 CHANNELS | 20 PROGRAMS

PROGRAM LIST

TENS AND PAIN

Conventional antalgic tens

Endorphinic tens

Rotator cuff tendinitis

Muscle pain

Knee osteoarthritis

post-surgical pain

Muscle injuries

Cervical pain

Trapezius pain

Sciatica

Lumbago

Epicondylitis

Scapulohumeral periarthritis

Carpal tunnel

Osteoarthritis

Menstrual pain

SPORT

Decontracting

Maximum strength

Cool down

FITNESS, BEAUTY, WELLNESS

Sculpting