

PREMIUM 400

4 CHANNELS ELECTROTHERAPY | 258 PROGRAMS

OVER 250 PROGRAMS IN YOUR HANDS, PREMIUM 400 IS THE IDEAL ELECTROSTIMULATOR FOR THE SPORT AND THE CARE OF THE BEAUTY, FOR MEN AND WOMEN.

58 FITNESS | 53 SPORT | 44 SPECIALS SPORT | 36 BEAUTY | 23 MICROCURRENTS | 18 SERIAL SEQUENTIAL STIMULATION
12 TENS AND PAIN | 7 ACTION NOW | 3 INCONTINENCE | 1 IONOPHORESIS

TECHNICAL FEATURES

Display

Visible area size 2,6"

Channels

4 independents (8 electrodes)

Frequency

0,3-150 Hz *

Pulse amplitude

40-450 µs *

Power

0-120 mA per channel

Power supply

Rechargeable batteries

Size

100x160x35 mm


Weight

450 gr

* According to the type of current

EQUIPMENT

- 1 Bag
- 1 Stimulator PREMIUM 400
- 4 Cables for electrodes connection
- 2 Cables for microcurrents and ionophoresis
- 4 Self-adhesive square electrodes
- 4 Self-adhesive rectangular electrodes
- 1 Charger
- 1 Operating manual



CE
0476



TENS



MENS



IONO



3S



ACTION NOW



RUN TIME



LAST 10



FAVORITES



MULTI USER



PROGRAMMABLE



STEM LOCK



FUNCTION 2+2



MY TRAINER



SYNCHRO STIM



WORK TIME



RECHARGEABLE BATTERIES

PREMIUM 400

4 CHANNELS ELECTROTHERAPY | 258 PROGRAMS

PROGRAMS LIST

SPORT

Capillarization

Warm-up

Pre-competition warm-up

Active recovery

Maximum strength

Endurance strength

Explosive strength

Aerobic endurance

Reactivity

Post-competition recovery

Decontracting

Hypertrophy

SPECIAL SPORT

Running

Cross-country skiing

Martial arts

Tennis

Soccer

Cycling

FITNESS

Firming

Bio-Pulse firming

Sculpting

Bio-Pulse sculpting

Toning

Mass building

Body sculpting

Definition

Jogging

Anaerobic fitness

Aerobic fitness

Cramp prevention

BEAUTY

Drainage

Bio-Pulse drainage

Lipolysis

Post-pregnancy lipolysis

Toning massage

Connective massage

Swollen arms

Face capillaries

Skin tone improvement

Post-pregnancy drainage

Post-pregnancy firming

Breast firming

Breast sculpting

Lifting effect

REHAB

Quadriceps atrophy

Recovery after ACL surgery

Shoulder subluxation prevention

TENS and PAIN

Conventional antalgic tens

Modulated antalgic tens

Endorphinic tens

Menstrual pain

Chronic pain

Cervical pain

Muscle pain

Knee pain

Scapulohumeral periarthritis

Chronic lumbago

Bursitis-tendinitis

Osteoarthritis

MICROCURRENTS

Epicondylitis

Periarthritis

Muscle restoration

Contusion

Hedema

Skin ulcers

Sciatica

Lumbago

Brachial neuralgia

Acute pain

Articular pain

Stiff neck

Whiplash

Cervical spondylosis

Shoulder sprain

Knee sprain

Carpal tunnel

Osteoarthritis

Ankle Sprain

Achille tendon inflammation

Patella tendon inflammation

Rotator cuff inflammation

Tendon Inflammation

INCONTINENCE

Mixed incontinence

Stress incontinence

URGE incontinence

IONOPHORESIS

Ionophoresis

SERIAL SEQUENTIAL STIMULATION

The "3S" program list includes 18 parameter combinations.

The "3S" programs are characterized by a delay activation of the channels 3 and 4 compared with the channels 1 and 2. The Serial Sequential Stimulation permits to stimulate the musculature in kinetic chain thanks to the differentiated activation times of the muscular groups involved.

ACTION NOW

The Action Now program list includes 7 parameter combinations.

This mode is suggested in sport field for athletic preparation where you wish to add the muscular contraction induced by a stimulator to a work made with overloads in dynamic and isometric form.