


KIT MAGNETOVET 4000 HORSE PRO




MAGNETOTHERAPY FOR HORSES

THE NEW MAGNETOVET 4000 HORSE PRO KIT, IN ADDITION TO THE DEVICE FOR MAGNETOTHERAPY WITH THE NEW PROGRAMS DEDICATED TO THE HORSE, INCLUDES ACCESSORIES DESIGNED FOR ITS TREATMENT: THERAPEUTIC BLANKET AND NECK-COVER, FLEXIBLE SOLENOIDS WITH RELATED SPLITTER CABLES, A PRACTICAL DEVICE CASE AND AN ADJUSTABLE STRAP FOR ATTACHING THE DEVICE TO THE HORSE. A CONVENIENT CARRYING BAG AND A DOUBLE BATTERY PACK COMPLETE THE KIT. MAGNETOVET 4000 HORSE IS THE HIGH INTENSITY PULSED ELECTROMAGNETIC FIELDS MAGNETOTHERAPY DEVICE SPECIALLY DESIGNED FOR THE TREATMENT OF HORSES. THE FOUR CHANNELS ALLOW THE SOLENOIDS, WHICH ARE INSERTED INTO THE THERAPEUTIC BLANKET AND NECK-COVER, TO BE CONNECTED. THIS ALLOWS THE ENTIRE SPINE OF THE HORSE TO BE TREATED.

THERAPEUTIC BLANKET




THERAPEUTIC NECK-COVER



STANDARD EQUIPEMENT | Code G6593

- 1 Therapeutic blanket
- 1 Therapeutic neck-cover
- 1 Case with shoulder strap for MAGNETOVET 4000 HORSE
- 1 Anti-shock carrying bag
- 1 Device MAGNETOVET 4000 HORSE
- 1 Battery pack (already included inside the device)
- 1 Magnum Test
- 1 Power supply unit
- 1 Operating manual
- 8 Flexible solenoids 40x15 cm for blanket and neck-cover
- 2 Splitting cables
- 1 Strap for protective case
- 1 Kit carry bag



ON REQUEST A WIDE RANGE OF ACCESSORIES COMPATIBLE WITH MAGNETOVET 4000 HORSE, IN ADDITION TO THOSE INSIDE THE KIT, IS AVAILABLE

KIT MAGNETOVET 4000 HORSE PRO

MAGNETOTHERAPY FOR HORSES

Rev. En_04_2024

DEVICE MAGNETOVET 4000 HORSE TECHNICAL FEATURES

Display

Size 2.6", backlit

Channels

4

Intensities

200 Gauss on average per channel, with flexible solenoids
800 Gauss on average as a whole (on 4 channels)

The nominal powers refer to accessories with two coils.
In the case of accessories with multiple coils, the power is evenly distributed

Frequency

Adjustable in programming from 5 to 200Hz

Duty Cycle

5% a 50%. Increment by steps of 5

Power supply

- Main
- Rechargeable batteries 4000mAh included

Size

17x22x6 cm

Weight

1,37 Kg

LIST OF PROGRAMS FOR HORSE

THE PROGRAMS ARE DIVIDED BY THE ANATOMICAL AREAS OF THE NECK, BACK AND LIMBS

Arthritis

Arthritis

Fracture

Delayed fracture consolidation

Discopathy

Tendinopathy

Wound

Ulcer

Anti-inflammatory

Injury

Contusion

Edema

Muscular relaxation

Decontracting

Muscular warm-up

Post-competition recovery

Muscle injury

LIST OF PROGRAMS FOR DOG/CAT

THE PROGRAMS INVOLVE THE USE OF THE THERAPEUTIC BED OR FLEXIBLE SOLENOIDS WITH 2 COILS, 30x10 cm

Fracture

Osteoarthritis

Discopathy

Cicatrization

Inflammatory pain

CE



- 2+2 MODE
- LAST 10
- PROGRAMMABLE
- BATTERY PACK
- MAINS POWER SUPPLY

KIT MAGNETOVET 4000 HORSE PRO

MAGNETOTHERAPY FOR HORSES

THERAPEUTIC BLANKET

TECHNICAL FEATURES

STRONG AND WASHABLE

Dimensions

142 cm lenght from withers (A) to tail (B)

Wieight

1,80 Kg (without solenoids)

Size

Universal

Belt

5 of adjustable fastening
3 of connecting belt to the therapeutic neck-cover

Pockets

8 for flexible solenoid insertion 40x15 cm, 4 coils
2 central for cable storage

Handles

3 for blanket positioning

Cables and extensions

Included

Material

96% polyester, 4% polypropylene



THERAPEUTIC NECK-COVER

TECHNICAL FEATURES

STRONG AND WASHABLE

Dimensions

121 cm length from ears (A) to withers (B)

Wieight

560 gr (without solenoids)

Size

Universal

Belts

3 fasteners, adjustable
3 connecting to the therapy blanket

Pockets

2 for inserting theflexible solenoids 40x15 cm with 4 coils

Cables and extensions

Included

Material

94% polyester, 3% polypropylene,
3% expanded polyethylene



KIT MAGNETOVET 4000 HORSE PRO

MAGNETOTHERAPY FOR HORSES

Rev. En_04_2024



PURCHASE OPTIONS

Cod. G6256 THERAPEUTIC FRONT BOOT

Cod. G6591 THERAPEUTIC FRONT BOOTS WITH PADDING AND SOLENOIDS, 2 COILS INCLUDED




Cod. G6255 THERAPEUTIC REAR BOOT

Cod. G6592 THERAPEUTIC REAR BOOTS WITH PADDING AND SOLENOIDS, 4 COILS INCLUDED



 **DOMINO S.R.L.**
Via Vittorio Veneto, 52
31013 Codognè (TV) - ITALY

 (+39) 0438 7933

 globuscorporation.com