

SonoCover

Ultrasound Probe Cover

The ultrasound covers, SonoCover, are either available in the linear standard version or with a pre-formed ultrasound probe chamber. The covers can be purchased optionally with a sterile 20ml ultrasound gel pack. Each package includes rubber bands and adhesive stripes for the fixation of the ultrasound probe cover.



Lateral Welded Seams

→ Artefact-free and accurate ultrasound image



Optional Ultrasound Probe Chamber

→ Pre-formed ultrasound probe chamber
→ Unhindered transmission of ultrasound waves



Telescope Folding

→ Simplifies the application of the cover under sterile conditions
→ Allows direct application of contact medium



Ultrasound Gel

→ Allows for best transmission of ultrasound waves

SonoCover

Ultrasound Probe Cover

SonoCover

with one stripe and two rubbers,
sterile

Size	Item no.	PU	Pre-formed chamber	Ultrasound gel 20ml
15 x 30 cm (6" x 12")	001151-1530	25		
15 x 60 cm (6" x 23 1/2")	001151-1560	25		
15 x 100 cm (6" x 40")	001151-1501	25		
15 x 30 cm (6" x 12")	011151-1530	25	●	
15 x 60 cm (6" x 23 1/2")	011151-1560	25	●	
15 x 100 cm (6" x 40")	011151-1501	25	●	
15 x 30 cm (6" x 12")	021151-1530	25		●
15 x 60 cm (6" x 23 1/2")	021151-1560	25		●
15 x 100 cm (6" x 40")	021151-1501	25		●
15 x 30 cm (6" x 12")	031151-1530	25	●	●
15 x 60 cm (6" x 23 1/2")	031151-1560	25	●	●
15 x 100 cm (6" x 40")	031151-1501	25	●	●

DISCLAIMER: Not all products are registered and approved for sale in all countries or regions. Indications of use may also vary by country and region. Please contact your country representative for product availability and information. Product images are for reference only.

PAJUNK® GmbH Medizintechnologie
Global Headquarters
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 9291-0
Fax: +49 7704 9291-600
info@pajunk.com
pajunk.com

PAJUNK® Medical Produkte GmbH
Central Europe
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 8008-0
Fax: +49 7704 8008-150
info@pajunk-medipro.com
pajunk.de/ pajunk.eu

PAJUNK® Swiss AG
Switzerland
Neuhofstraße 5A
6340 Baar, Switzerland
Tel.: +41 41 760 96 27
info@pajunk.ch
pajunk.ch

PAJUNK® Medical System L.P.
United States of America
4575 Marconi Drive
Alpharetta, GA 30005, USA
Tel.: +1 888 972 5865
Fax: +1 678 514 3388
info@pajunk-usa.com
pajunkusa.com

PAJUNK® UK Medical Products Ltd
United Kingdom
Unit C1 The Waterfront
Goldcrest Way, Newburn Riverside
Newcastle upon Tyne, NE15 8NY, UK
Tel.: +44 191 264 7333
info@pajunk.co.uk
pajunk.co.uk