



Meril

ENT



Breathe freely

MeticTM

Airway Dilation Catheter

The Metic™ Balloon Dilation System is designed to safely dilate airway strictures with minimal mucosal trauma. Developed in partnership with ENT Surgeons, the Metic™ offers an safe, elegant, minimally invasive solution for airway management.

Patient safety is our commitment.

The Metic™ delivers:

A product specifically designed for the airway

- Custom balloon catheter and stylet optimized for airway anatomy

Atraumatic access

- Low-profile system crosses narrow stenoses
- Soft seeker tip guides device through delicate anatomy


Precise placement under direct visualization

- Integrated stylet provides enhanced support and maneuverability
- Malleable shaft can be custom-shaped for easy access

Controlled radial dilation

- Durable high pressure balloon engineered for resistant lesions
- Non-compliant design maintains consistent diameter through multiple dilations
- Precise inflation and pressure regulation

Airway Balloon Dilation




Atraumatic access

Safe, accurate placement



controlled radial dilation

Consistent dilation of strictures



Catalog Number	Balloon Diameter (mm)	Balloon Length (mm)	Units (per box)
BC0524A	5	24	5
BC0724A	7	24	5
BC0924A	9	24	5
BC1040A	10	40	5
BC1240A	12	40	5
BC1440A	14	40	5
BC1640A	16	40	5