

e-xam

Diagnostic penlight



Photobiologically validated high-performance pen-lights for medical diagnosis of throat and pupil. Available with LED or Xenon illumination.

- Photobiologically validated for lighting of pupils according to EN 62471
- Pre-focused and even field of light
- Bright and natural white light by high-performance LED
- LED burns at least 20,000 hours
- Low maintenance costs thanks to long LED and battery service lives
- ABS casing with metal clip, disinfectable easy to use with ON/OFF switch including blade holder, batteries and nylon case
- Designed for durability, comes with a 2-year warranty including LED



Technical data:

e-xam LED 2.5 V	5000 - 5500 K 180,000 Lux max. field of light: approx. 11 cm diameter (distance 15 cm)
e-xam Xenon 2.5 V	3500 - 4000 K 60,000 Lux max. field of light: approx. 14 cm diameter (distance 15 cm)
Batteries	2 x Alkaline AA, 1.5 V
5130-01	e-xam, XL 2.5 V, black, with tongue blade holder
5131-01	e-xam, LED 2.5 V, black, with tongue blade holder

Please visit **riester.de** for more information