

ri-pen®

Diagnostic pupil penlight



Sustainable penlight for medical diagnosis with LED illumination

- ▶ Photobiologically validated for lighting of pupils according to EN 62471
- ▶ New design with pupil gauge from 2 mm - 9 mm and a ruler with 8 cm [3 in]
- ▶ UV free light through LED illumination in an aluminium casing with practical metal clip
- ▶ Supplied as single units or in packs of 6 (assorted colors or 6 per color)
- ▶ Simple ON/OFF function with push button on the end cap of the penlight
- ▶ Comes with 2 AAA batteries at time of purchase

Technical data:

ri-pen® LED	3 V		
	5300 K		
	40,000 Lux max.		
Field of view	Approximately 40 mm diameter at a working distance of 15 cm		
Batteries	2 x Alkaline AAA, 1.5 V		
5075-526	ri-pen® black, pack of 1 pc.	5075	ri-pen® black, pack of 6 pcs.
5074-526	ri-pen® silver, pack of 1 pc.	5074	ri-pen® silver, pack of 6 pcs.
5077-526	ri-pen® red, pack of 1 pc.	5077	ri-pen® red, pack of 6 pcs.
5078-526	ri-pen® purple, pack of 1 pc.	5078	ri-pen® purple, pack of 6 pcs.
5071-526	ri-pen® blue, pack of 1 pc.	5071	ri-pen® blue, pack of 6 pcs.
5072-526	ri-pen® green, pack of 1 pc.	5072	ri-pen® green, pack of 6 pcs.
		5070	ri-pen® assorted colors, pack of 6 pcs.

Please visit **riester.de** for more information