

gymna[®]

TECHNICAL DATASHEET

ShockMaster 300



Technical device specifications

Unit dimensions (w x d x h)	36 x 38 x 16 cm
Unit weight (kg)	11,9 kg
Connection channels	1 channel
Touchscreen	7 inch Full Colour, 180° view
Frequency range	0.5 - 17 Hz (in steps of 0.1 Hz)
Maximum frequency range with V-Actor	31 Hz
Pressure range	0,3 - 4 bar _{eff} (in steps of 0.1 bar _{eff})
Classification according to MDD 93/42/EEC	Class IIa
Classification according to MDR 509-281	Class II
Electrical protection class	Class I
Applied parts type	Type B
Ingress protection	IPX1
Notified body	CE0344

Technical handpiece specifications

On / Off function	Electronic start / stop button No air cooling
Weight (g)	+/- 500 g
Ergonomics	Integrated soft rubber protection Ergonomic shape Shockabsorbing rubber grip parts



device



handpiece

Available marketing material

Waiting room poster certificate

An eye-catching and professional poster designed to inform and reassure patients while they wait, highlighting the benefits of ShockMaster therapy and your ShockMaster expertise.

Application guide (digital)

A clear, step-by-step digital guide with treatment protocols and practical tips for effective daily use.

Patient folder (digital)

A patient-friendly digital folder providing accessible information about treatments, benefits, and what to expect.

Everything you
need to showcase
the ShockMaster
with confidence.

Ordering information

Gymna ShockMaster 300

342.300	Gymna ShockMaster 300
---------	-----------------------

Standard included: Device Gymna ShockMaster 300, RSW handpiece Precision*, beam applicator, deep impact applicator, D-actor 20mm applicator, 4 revision sets (2 projectiles and tubes inside handpiece box + 1 revision kit handpiece), 1 protective antiskid mat, 1 bottle Physio Care contact neutral hydrogel 500 ml, 1 tissuebox, 1 waiting room poster, 1 certificate and 1 application guide, ring binder with documents and information.

Gymna ShockMaster Accessories

342.720	Gymna ShockMaster RSW handpiece Precision
342.716	Gymna ShockMaster RSW handpiece short
342.749	Gymna ShockMaster body replacement set RSW handpiece short
309.771	Gymna ShockMaster revision kit RSW handpiece
342.705	Gymna ShockMaster revision kit RSW handpiece short / handpiece Precision
338.239	Gymna ShockMaster V-actor
342.727	Gymna ShockMaster fascia applicator set
342.738	Gymna ShockMaster spine applicator set
338.613	Gymna ShockMaster classic 15 mm applicator

*In official documentation, the handpiece Precision is referred to as: ShockMaster RSW handpiece Short W

Ordering information

338.618	Gymna ShockMaster beam applicator
338.624	Gymna ShockMaster deep impact applicator
338.630	Gymna ShockMaster focus lens applicator
342.760	Gymna ShockMaster atlas applicator
338.602	Gymna ShockMaster acupuncture 6 mm applicator
342.771	Gymna ShockMaster ceramic applicator
342.782	Gymna ShockMaster trigger 10 mm applicator
338.646	Gymna ShockMaster D-actor 20 mm applicator
338.657	Gymna ShockMaster D-actor 35 mm applicator
342.793	Gymna ShockMaster D-actor 20 mm "T" applicator
338.668	Gymna ShockMaster V-actor 25 mm applicator
338.679	Gymna ShockMaster V-actor 40 mm applicator
338.217	Gymna ShockMaster protective anti skid mat
342.551	Gymna ShockMaster poster
320.804	Gymna Mobile equipment trolley
338.316	Gymna ShockMaster custom made carrying bag
342.673	Gymna ShockMaster application guide EN
341.088	Gymna Physio Care Contact neutral hydrogel 500 ml
341.099	Gymna Physio Care Contact neutral hydrogel 5 l

